

ZMPSTE24 Positive Control for STJ503602 peptide (STJ505266)

STJ505266

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ZMPSTE24 Positive Control for STJ503602 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10269
Gene Symbol	ZMPSTE24
Uniprot ID	FACE1_HUMAN
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with ZMPSTE24 antibody STJ503602.
Immunogen Sequence	