

PDE4B Blocking Peptide for STJ502320 peptide (STJ505276)

STJ505276

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE4B Blocking Peptide for STJ502320 is synthetically produced from the C-Term sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5142
Gene Symbol	PDE4B
Uniprot ID	PDE4B_HUMAN
Immunogen	Synthetic peptide common to all PDE4B subtypes taken near C-terminal
Immunogen Region	C-Term
Specificity	This blocking peptide is recommended for use in combination with PDE4B antibody, STJ502320
Immunogen Sequence	