

PDE11A Blocking Peptide peptide (Phospho) (S162) (STJ505298)

STJ505298

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PDE11A Blocking Peptide is synthetically produced from the 100-200 sequence about the S162 and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	50940
Gene Symbol	PDE11A
Uniprot ID	PDE11_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-200 on PDE11A protein. Phosphorylated Serine 162 residue.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	