

PDE12 Blocking Peptide for STJ502225 peptide (STJ505305)

STJ505305

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE12 Blocking Peptide for STJ502225 is synthetically produced from the 570-620 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	201626
Gene Symbol	PDE12
Uniprot ID	PDE12_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 570-620 on human PDE12A protein.
Immunogen Region	570-620
Specificity	This blocking peptide is recommended for use in combination with PDE12 antibody, STJ502225
Immunogen Sequence	