

ZDHHC9 Blocking Peptide for STJ503571 peptide (STJ505313)

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ZDHHC9 Blocking Peptide for STJ503571 is synthetically produced from the 250-300 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51114
Gene Symbol	ZDHHC9
Uniprot ID	ZDHC9_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300 on human Palmitoyltransferase ZDHHC9 protein.
Immunogen Region	250-300
Specificity	This blocking peptide is recommended for use in combination with ZDHHC9 antibody, STJ503571
Immunogen Sequence	