

Humanin Positive Control for STJ506139 peptide (STJ505327)

STJ505327

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Humanin Positive Control for STJ506139 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	WB: 1:500
Range	ELISA: 1:10, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	MT-RNR2
Uniprot ID	HUNIN_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	This is positive control is recommended for use in combination with Humanin antibody STJ506139.
Immunogen	
Sequence	