

PDE2A2 Blocking Peptide peptide (STJ505340)

STJ505340

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE2A2 Blocking Peptide is synthetically produced from the 1-50 aa sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5138
Gene Symbol	PDE2A
Uniprot ID	PDE2A_HUMAN
Immunogen	Synthetic peptide corresponding to region within amino acids 1-50 on N-terminus of human PDE2A2 protein.
Immunogen Region	1-50 aa
Specificity	
Immunogen Sequence	