

MeCP2 isoform 1 Positive Control for STJ501740 peptide (STJ505350)

STJ505350

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MeCP2 isoform 1 Positive Control for STJ501740 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4204
Gene Symbol	MECP2
Uniprot ID	MECP2_HUMAN
Immunogen	Purified Protein
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with MeCP2 isoform 1 antibody STJ501740.
Immunogen Sequence	