

PDE7A Blocking Peptide peptide (Phospho) (STJ505351)

STJ505351

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PDE7A Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	18583
Gene Symbol	Pde7a
Uniprot ID	PDE7A_MOUSE
Immunogen	Synthetic cyclic peptide common to PDE7A1 and PDE7A2 variants.
Immunogen Region	
Specificity	
Immunogen	
Sequence	