

PDE7A Positive Control peptide (Phospho) (STJ505352)

STJ505352

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PDE7A Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	18583
Gene Symbol	Pde7a
Uniprot ID	PDE7A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	