

FGFR1 Blocking Peptide peptide (Phospho) (STJ505359)

STJ505359

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-FGFR1 Blocking Peptide is synthetically produced from the 640-690 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2260
Gene Symbol	FGFR1
Uniprot ID	FGFR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 640-690 on human basic fibroblast growth factor receptor protein 1.
Immunogen Region	640-690
Specificity	
Immunogen	
Sequence	