

PAX2 isoform a + e Blocking Peptide for STJ502140 peptide (STJ505371)

STJ505371

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PAX2 isoform a + e Blocking Peptide for STJ502140 is synthetically produced from the 200-220 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5076
Gene Symbol	PAX2
Uniprot ID	PAX2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-220 on human PAX2 isoform d protein.
Immunogen Region	200-220
Specificity	This blocking peptide is recommended for use in combination with PAX2 isoform a + e antibody, STJ502140
Immunogen Sequence	