

OGT1 Blocking Peptide peptide (STJ505421)

STJ505421

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	OGT1 Blocking Peptide is synthetically produced from the 980-1030 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8473
Gene Symbol	OGT
Uniprot ID	OGT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 980-1030 on human O-linked N-acetylglucosamine transferase protein.
Immunogen Region	980-1030
Specificity	
Immunogen Sequence	