

MAPK4 Blocking Peptide for STJ501706 peptide (STJ505444)

STJ505444

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MAPK4 Blocking Peptide for STJ501706 is synthetically produced from the 500-550 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5596
Gene Symbol	MAPK4
Uniprot ID	MK04_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human Mitogen-activated protein kinase 4.
Immunogen Region	500-550
Specificity	This blocking peptide is recommended for use in combination with MAPK4 antibody, STJ501706
Immunogen Sequence	