

ABHD2 Positive Control for STJ500017 peptide (STJ505473)

STJ505473

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Biomolecules |
| Short Description | ABHD2 Positive Control for STJ500017 is synthetically produced from the sequence and is suitable for use in western blot applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: 1:500 |
| Formulation | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 54608 |
| Gene Symbol | Abhd2 |
| Uniprot ID | ABHD2_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | This is positive control is recommended for use in combination with ABHD2 antibody STJ500017. |
| Immunogen Sequence | |