

5HT2C Receptor Blocking Peptide peptide (STJ505510)

STJ505510

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	5HT2C Receptor Blocking Peptide is synthetically produced from the 400-458 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3358
Gene Symbol	HTR2C
Uniprot ID	5HT2C_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-458 on human 5-hydroxytryptamine receptor 2C.
Immunogen Region	400-458
Specificity	
Immunogen Sequence	