

## IKAP Positive Control for STJ501408 peptide (STJ505519)

STJ505519

### GENERAL INFORMATION

|                          |                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Biomolecules                                                                                                                          |
| <b>Short Description</b> | IKAP Positive Control for STJ501408 is synthetically produced from the sequence and is suitable for use in western blot applications. |
| <b>Applications</b>      | WB                                                                                                                                    |
| <b>Host/Source</b>       |                                                                                                                                       |
| <b>Reactivity</b>        |                                                                                                                                       |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      |                                                                                                                                                                                                                                                          |
| <b>Clone ID</b>       |                                                                                                                                                                                                                                                          |
| <b>Concentration</b>  |                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>    |                                                                                                                                                                                                                                                          |
| <b>Purification</b>   |                                                                                                                                                                                                                                                          |
| <b>Dilution Range</b> | WB: 1:500                                                                                                                                                                                                                                                |
| <b>Formulation</b>    | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| <b>Isotype</b>        |                                                                                                                                                                                                                                                          |
| <b>Storage</b>        | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                                                                                                          |
| <b>Instruction</b>    |                                                                                                                                                                                                                                                          |

### TARGET INFORMATION

|                           |                                                                                              |
|---------------------------|----------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">230233</a>                                                                       |
| <b>Gene Symbol</b>        | <a href="#">Elp1</a>                                                                         |
| <b>Uniprot ID</b>         | <a href="#">ELP1_MOUSE</a>                                                                   |
| <b>Immunogen</b>          |                                                                                              |
| <b>Immunogen Region</b>   |                                                                                              |
| <b>Specificity</b>        | This is positive control is recommended for use in combination with IKAP antibody STJ501408. |
| <b>Immunogen Sequence</b> |                                                                                              |