

## Syncytin 1 Blocking Peptide peptide (STJ505530)

STJ505530

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Syncytin 1 Blocking Peptide is synthetically produced from the 480-538 sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">30816</a>
<b>Gene Symbol</b>	<a href="#">ERVW-1</a>
<b>Uniprot ID</b>	<a href="#">SYCY1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 480-538 on human Syncytin-1 protein.
<b>Immunogen Region</b>	480-538
<b>Specificity</b>	
<b>Immunogen Sequence</b>	