

Syncytin 2 Blocking Peptide peptide (STJ505533)

STJ505533

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Syncytin 2 Blocking Peptide is synthetically produced from the 80-130 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	405754
Gene Symbol	ERVFRD-1
Uniprot ID	SYCY2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 80-130 on human Syncytin-2 protein.
Immunogen Region	80-130
Specificity	
Immunogen Sequence	