

IKAP Blocking Peptide for STJ501408 peptide (STJ505537)

STJ505537

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	IKAP Blocking Peptide for STJ501408 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	230233
Gene Symbol	Elp1
Uniprot ID	ELP1_MOUSE
Immunogen	Recombinant IkappaB kinase complex-associated protein.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with IKAP antibody, STJ501408
Immunogen Sequence	