

D4S234E Blocking Peptide for STJ501974 peptide (STJ505538)

STJ505538

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	D4S234E Blocking Peptide for STJ501974 is synthetically produced from the 140-190 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25247
Gene Symbol	Nsg1
Uniprot ID	NSG1_RAT
Immunogen	Synthetic peptide taken within amino acid region 140-190 on human Neuron-specific protein family member 1.
Immunogen Region	140-190
Specificity	This blocking peptide is recommended for use in combination with D4S234E antibody, STJ501974
Immunogen Sequence	