

D4S234E Positive Control for STJ501973 peptide (STJ505539)

STJ505539

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	D4S234E Positive Control for STJ501973 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25247
Gene Symbol	Nsg1
Uniprot ID	NSG1_RAT
Immunogen	Purified Protein
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with D4S234E antibody STJ501973.
Immunogen Sequence	