

Von Hippel Lindau Positive Control peptide (STJ505541)

STJ505541

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Von Hippel Lindau Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7428
Gene Symbol	VHL
Uniprot ID	VHL_HUMAN
Immunogen	Purified Protein
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	