

NSG2 Blocking Peptide for STJ501977 peptide (STJ505542)

STJ505542

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	NSG2 Blocking Peptide for STJ501977 is synthetically produced from the 130-171 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51617
Gene Symbol	NSG2
Uniprot ID	NSG2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 130-171 on human Neuron-specific protein family member 2.
Immunogen Region	130-171
Specificity	This blocking peptide is recommended for use in combination with NSG2 antibody, STJ501977
Immunogen Sequence	