

Deformed Epidermal Autoregulatory Factor 1 Blocking Peptide for STJ500748 peptide (STJ505547)

STJ505547

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Deformed Epidermal Autoregulatory Factor 1 Blocking Peptide for STJ500748 is synthetically produced from the 500-565 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	54006
Gene Symbol	Deaf1
Uniprot ID	DEAF1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 500-565 on mouse Deformed epidermal autoregulatory factor 1 homolog protein.
Immunogen Region	500-565
Specificity	This blocking peptide is recommended for use in combination with Deformed Epidermal Autoregulatory Factor 1 antibody, STJ500748
Immunogen Sequence	