

THAP1 Positive Control for STJ503247 peptide (STJ505550)

STJ505550

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	THAP1 Positive Control for STJ503247 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6548
Gene Symbol	SLC9A1
Uniprot ID	SL9A1_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with THAP1 antibody STJ503247.
Immunogen Sequence	