

TSG101 Blocking Peptide for STJ503421 peptide (STJ505553)

STJ505553

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TSG101 Blocking Peptide for STJ503421 is synthetically produced from the 330-390 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7251
Gene Symbol	TSG101
Uniprot ID	TS101_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 330-390 on Tumor susceptibility gene 101 protein.
Immunogen Region	330-390
Specificity	This blocking peptide is recommended for use in combination with TSG101 antibody, STJ503421
Immunogen Sequence	