

Cleavage stimulation factor subunit 1 Positive Control peptide (STJ505580)

STJ505580

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Cleavage stimulation factor subunit 1 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1477
Gene Symbol	CSTF1
Uniprot ID	CSTF1_HUMAN
Immunogen	Purified Protein
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	