

HIPK3 Positive Control for STJ501347 peptide (STJ505592)

STJ505592

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	HIPK3 Positive Control for STJ501347 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10114
Gene Symbol	HIPK3
Uniprot ID	HIPK3_HUMAN
Immunogen	Purified Protein
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with HIPK3 antibody STJ501347.
Immunogen Sequence	