

AMPK alpha 1 Positive Control for STJ500075 peptide (STJ505611)

STJ505611

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	AMPK alpha 1 Positive Control for STJ500075 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65248
Gene Symbol	Prkaa1
Uniprot ID	AAPK1_RAT
Immunogen	Purified Protein
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with AMPK alpha 1 antibody STJ500075.
Immunogen Sequence	