

MeCP2 isoform 1 Blocking Peptide for STJ501738 peptide (STJ505640)

STJ505640

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MeCP2 isoform 1 Blocking Peptide for STJ501738 is synthetically produced from the 1-50 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	17257
Gene Symbol	Mecp2
Uniprot ID	MECP2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1-50 on mouse Methyl-CpG-binding Protein 2 isoform 1 protein.
Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with MeCP2 isoform 1 antibody, STJ501738
Immunogen Sequence	