

Twist Blocking Peptide peptide (STJ505655)

STJ505655

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Twist Blocking Peptide is synthetically produced from the 50-100 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7291
Gene Symbol	TWIST1
Uniprot ID	TWIST1_HUMAN
Immunogen	Synthetic peptide take within amino acid region 50-100 on human Twist-related protein 1.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	