

Elf4/MEF Blocking Peptide for STJ500869 peptide (STJ505657)

STJ505657

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Elf4/MEF Blocking Peptide for STJ500869 is synthetically produced from the 605-655 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56501
Gene Symbol	Elf4
Uniprot ID	ELF4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 605-655 on mouse ETS-related transcription factor Elf-4 protein.
Immunogen Region	605-655
Specificity	This blocking peptide is recommended for use in combination with Elf4/MEF antibody, STJ500869
Immunogen Sequence	