

## Estrogen Related Receptor gamma Positive Control peptide (STJ505660)

STJ505660

### GENERAL INFORMATION

|                          |                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Biomolecules                                                                                                                                       |
| <b>Short Description</b> | Estrogen Related Receptor gamma Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications. |
| <b>Applications</b>      | WB                                                                                                                                                 |
| <b>Host/Source</b>       |                                                                                                                                                    |
| <b>Reactivity</b>        |                                                                                                                                                    |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      |                                                                                                                                                                                                                                                           |
| <b>Clone ID</b>       |                                                                                                                                                                                                                                                           |
| <b>Concentration</b>  |                                                                                                                                                                                                                                                           |
| <b>Conjugation</b>    |                                                                                                                                                                                                                                                           |
| <b>Purification</b>   |                                                                                                                                                                                                                                                           |
| <b>Dilution Range</b> | WB: 1:2, 500-1:5, 000                                                                                                                                                                                                                                     |
| <b>Formulation</b>    | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| <b>Isotype</b>        |                                                                                                                                                                                                                                                           |
| <b>Storage</b>        | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                                                                                                           |
| <b>Instruction</b>    |                                                                                                                                                                                                                                                           |

### TARGET INFORMATION

|                    |                                                   |
|--------------------|---------------------------------------------------|
| <b>Gene ID</b>     | 2104                                              |
| <b>Gene Symbol</b> | ESRRG                                             |
| <b>Uniprot ID</b>  | ERR3_HUMAN                                        |
| <b>Immunogen</b>   | Purified Estrogen Related receptor gamma protein. |
| <b>Immunogen</b>   |                                                   |
| <b>Region</b>      |                                                   |
| <b>Specificity</b> |                                                   |
| <b>Immunogen</b>   |                                                   |
| <b>Sequence</b>    |                                                   |