

**Brk Positive Control for STJ502711 peptide (STJ505666)**

STJ505666

**GENERAL INFORMATION**

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Brk Positive Control for STJ502711 is synthetically produced from the sequence and is suitable for use in western blot applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:2, 500-1:5, 000
<b>Formulation</b>	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">20459</a>
<b>Gene Symbol</b>	<a href="#">Ptk6</a>
<b>Uniprot ID</b>	<a href="#">PTK6_MOUSE</a>
<b>Immunogen</b>	Purified Brk protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	This is positive control is recommended for use in combination with Brk antibody STJ502711.
<b>Immunogen Sequence</b>	