

## BLK Blocking Peptide peptide (STJ505671)

STJ505671

### GENERAL INFORMATION

|                          |                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Biomolecules                                                                                                           |
| <b>Short Description</b> | BLK Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications. |
| <b>Applications</b>      | Immunodepletion/Immunocompetition                                                                                      |
| <b>Host/Source</b>       |                                                                                                                        |
| <b>Reactivity</b>        |                                                                                                                        |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                                |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      |                                                                                                                                                                |
| <b>Clone ID</b>       |                                                                                                                                                                |
| <b>Concentration</b>  |                                                                                                                                                                |
| <b>Conjugation</b>    |                                                                                                                                                                |
| <b>Purification</b>   |                                                                                                                                                                |
| <b>Dilution Range</b> |                                                                                                                                                                |
| <b>Formulation</b>    | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| <b>Isotype</b>        |                                                                                                                                                                |
| <b>Storage</b>        | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                |
| <b>Instruction</b>    |                                                                                                                                                                |

### TARGET INFORMATION

|                           |                                     |
|---------------------------|-------------------------------------|
| <b>Gene ID</b>            | <a href="#">640</a>                 |
| <b>Gene Symbol</b>        | <a href="#">BLK</a>                 |
| <b>Uniprot ID</b>         | <a href="#">BLK_HUMAN</a>           |
| <b>Immunogen</b>          | Tyrosine-protein kinase BLK Protein |
| <b>Immunogen Region</b>   |                                     |
| <b>Specificity</b>        |                                     |
| <b>Immunogen Sequence</b> |                                     |