

RBP2 Blocking Peptide for STJ502757 peptide (STJ505683)

STJ505683

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	RBP2 Blocking Peptide for STJ502757 is synthetically produced from the 200-250 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	74152
Gene Symbol	Stra6l
Uniprot ID	STR6L_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse Retinol-Binding Protein 4 Receptor 2/RBP4 receptor 2/STRA6-Like Protein.
Immunogen Region	200-250
Specificity	This blocking peptide is recommended for use in combination with RBP2 antibody, STJ502757
Immunogen Sequence	