

TMC1 Positive Control for STJ503319 peptide (STJ505686)

STJ505686

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TMC1 Positive Control for STJ503319 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	117531
Gene Symbol	TMC1
Uniprot ID	TMC1_HUMAN
Immunogen	Purified Transmembrane channel-like protein 1.
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with TMC1 antibody STJ503319.
Immunogen Sequence	