

Filamin B Blocking Peptide peptide (Phospho) (2481) (STJ505693)

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Filamin B Blocking Peptide is synthetically produced from the 2450-2500 sequence about the 2481 and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2317
Gene Symbol	FLNB
Uniprot ID	FLNB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 2450-2500 on human Phospho-Filamin-B. Phosphorylated 2481 serine residue.
Immunogen	2450-2500
Region	
Specificity	
Immunogen	
Sequence	