

PTP alpha Positive Control peptide (STJ505707)

STJ505707

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PTP alpha Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5786
Gene Symbol	PTPRA
Uniprot ID	PTPRA_HUMAN
Immunogen	Receptor-type tyrosine-protein phosphatase Alpha
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	