

PICK1 Positive Control for STJ505933 peptide (STJ505721)

STJ505721

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PICK1 Positive Control for STJ505933 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9463
Gene Symbol	PICK1
Uniprot ID	PICK1_HUMAN
Immunogen	Purified PRKCA-binding protein/PICK1.
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with PICK1 antibody STJ505933.
Immunogen Sequence	