

PP2A-beta Blocking Peptide for STJ505947 peptide (STJ505730)

STJ505730

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PP2A-beta Blocking Peptide for STJ505947 is synthetically produced from the 450-497 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5526
Gene Symbol	PPP2R5B
Uniprot ID	2A5B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-497 on human Serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Beta subunit.
Immunogen Region	450-497
Specificity	This blocking peptide is recommended for use in combination with PP2A-beta antibody, STJ505947
Immunogen Sequence	