

PP2A-delta Blocking Peptide for STJ505950 peptide (STJ505732)

STJ505732

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PP2A-delta Blocking Peptide for STJ505950 is synthetically produced from the 540-600 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5528
Gene Symbol	PPP2R5D
Uniprot ID	2A5D_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 540-600 on human Serine/Threonine-Protein Phosphatase 2A B56-Delta subunit.
Immunogen Region	540-600
Specificity	This blocking peptide is recommended for use in combination with PP2A-delta antibody, STJ505950
Immunogen Sequence	