

PP2A-delta Positive Control for STJ505952 peptide (STJ505733)

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Biomolecules |
| Short Description | PP2A-delta Positive Control for STJ505952 is synthetically produced from the sequence and is suitable for use in western blot applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 5528 |
| Gene Symbol | PPP2R5D |
| Uniprot ID | 2A5D_HUMAN |
| Immunogen | Purified Serine/Threonine-Protein Phosphatase 2A B56-Delta subunit. |
| Immunogen Region | |
| Specificity | This is positive control is recommended for use in combination with PP2A-delta antibody STJ505952. |
| Immunogen Sequence | |