

PP2A Activator protein Blocking Peptide for STJ505961 peptide (STJ505738)

STJ505738

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PP2A Activator protein Blocking Peptide for STJ505961 is synthetically produced from the 200-250 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Serine/Threonine-Protein Phosphatase 2A Activator.
Immunogen Region	200-250
Specificity	This blocking peptide is recommended for use in combination with PP2A Activator protein antibody, STJ505961
Immunogen Sequence	