

Bruchpilot Blocking Peptide peptide (Phospho) (STJ505742)

STJ505742

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Bruchpilot Blocking Peptide is synthetically produced from the 150-200 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 150-200 on D. melanogaster (fruit fly) Bruchpilot protein.
Immunogen Region	150-200
Specificity	
Immunogen	
Sequence	