

CRALBP Blocking Peptide for STJ505978 peptide (STJ505750)

STJ505750

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	CRALBP Blocking Peptide for STJ505978 is synthetically produced from the 260-317 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6017
Gene Symbol	RLBP1
Uniprot ID	RLBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-317 on human Retinaldehyde-binding protein 1/CRALBP/RLBP1 protein.
Immunogen Region	260-317
Specificity	This blocking peptide is recommended for use in combination with CRALBP antibody, STJ505978
Immunogen Sequence	