

c-Fos Positive Control for STJ506017 peptide (STJ505777)

STJ505777

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	c-Fos Positive Control for STJ506017 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2353
Gene Symbol	FOS
Uniprot ID	FOS_HUMAN
Immunogen	Purified Proto-oncogene c-Fos.
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with c-Fos antibody STJ506017.
Immunogen Sequence	