

PBXIP1 Blocking Peptide peptide (Phospho) (STJ505778)

STJ505778

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PBXIP1 Blocking Peptide is synthetically produced from the 100-150 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	229534
Gene Symbol	Pbxip1
Uniprot ID	PBIP1_MOUSE
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 100-150 on mouse Pre-B-cell leukemia transcription factor-interacting protein 1.
Immunogen	100-150
Region	
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081