

Hck Blocking Peptide for STJ506037 peptide (STJ505788)

STJ505788

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Hck Blocking Peptide for STJ506037 is synthetically produced from the 350-400 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3055
Gene Symbol	HCK
Uniprot ID	HCK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human Hematopoietic cell kinase/Tyrosine-protein kinase HCK.
Immunogen Region	350-400
Specificity	This blocking peptide is recommended for use in combination with Hck antibody, STJ506037
Immunogen Sequence	